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BALANZAS DE PRECISIÓN

---

# AND-GR

## 200-202

**Moretti**

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Quick Start Guide – GR Series Analytical Balances .....	1
Guide de démarrage rapide – Série GR Balances analytiques .....	13
Guía de inicio rápido – Serie GR Básculas analíticas .....	25


# GR Series Analytical Balances

## Quick Start Guide



 This is a hazard alert mark

# UNPACKING

- Unpack the balance carefully and keep the packing material if you want to transport the balance.
- In the packing box you will find this manual and the following:

Step Cards (A step card is attached to the balance)

---

### SETTING UP

---



To ensure that you get the best performance from your balance, please observe the following conditions as closely as possible:

- The best operating temperature is about 20°C / 68°F at about 50% Relative Humidity.
- Try to ensure a stable power source when using the AC adapter.
- Please warm-up the balance for at least one hour. Plug-in the AC adapter as usual.
- The weighing room should be free of dust.
- The weighing table should be solid and free from vibration, drafts (from frequently opening doors or windows) and as level as possible.
- Keep the balance level by using the bubble spirit level.
- Do not install the balance near heaters or air conditioners.
- Do not install the balance in direct sunlight.
- Do not use the balance near other equipment which produces magnetic fields.
- Corners of rooms are best as they are less prone to vibrations.
- Calibrate the balance before operation and after moving it to another location.
- Do not place or use the balance where there is flammable or corrosive gas present.
- Do not use solvents to clean the balance. For best cleaning, wipe with a dry, lint free cloth or a lint free cloth that is moistened with warm water and a mild detergent.
- When you transport the balance, hold it as shown in the illustration. Never lift the balance using the weighing chamber frame.
- Keep magnetic substances away from the balance.
- Avoid mechanical shock to your balance.
- Avoid dust and water so that the balance weighs correctly. Protect the internal parts from liquid spills and excessive dust.
- Remove and clean the floor plate of the weighing chamber.
- Use the special shipping box supplied for transportation.
- Do not disassemble the balance. Contact your local A&D dealer if your balance needs service or repair.

---

## PART NAMES AND FUNCTIONS

---

# GR Series

## Analytical Balances

## Quick Start Guide

There are two operation types for pressing a key. Each key operation performs a different function:

First type: "Press and release the key immediately" or "Click the key"

Second type: "Press and hold the key"

The first type is "to press the key". This first type is normal key operation during measurement.

Do not press and hold the key if you do not want to perform a rewrite of the internal parameters.

Press the key

(Press and release the key immediately)

Press and hold the key

Key	Press the key	Press and hold the key
ON:OFF 	Display ON / OFF key. The standby indicator is displayed when the balance is turned off with this key. Weighing data is displayed when the balance is turned on with this key.	
RANGE 	Minimum figure of weighing data is changed.	The function table menu is displayed. Refer to section 9 of the Instruction Manual "Function Table"
MODE 	Units are changed (selected from the function table). Refer to section 4 of the Instruction Manual "Weighing Units".	Response adjustment is performed.
CAL 	This key performs calibration of the balance using the internal weight.	Other items of the calibration menu are displayed.
PRINT 	Weighing data is stored in the balance (factory setting) or is output to the RS-232C interface. This key functions according to the function table.	Data memory menu or GLP menu is displayed. This key functions according to the function table. Factory setting is "not used".
RE-ZERO 	The key sets the display to zero. This key returns a weighing value to the centre of zero when the weighing pan is empty and can also tare (cancel) the weight of container and/or sample. Please use this key before each weighing to cancel possible error.	

---

## CALIBRATION

---

### Calibration Using the Internal mass

1. Connect the AC adaptor and warm up the balance for at least one hour with nothing on the pan.
2. Press and hold the CAL key until the balance displays "cc in" and then release the key.
3. The balance measures the zero point. Prevent vibration and drafts from affecting the balance.
4. The measured zero-point data is displayed.
5. Ready for the internal weight measurement.
6. The balance measures the internal weight. Prevent vibration and drafts from affecting the balance.
7. The internal weight data is displayed.
8. The balance informs you when the calibration test is finished. If the "GLP output (info)" of the "Function Table" is set to "1" or "2", "Calibration Test Report" is output by the RS-232C interface.
9. The balance will automatically return to the weighing mode after the calibration test is finished.

### Calibration Using the External mass




1. Connect the AC adaptor and warm up the balance for at least one hour with nothing on the pan.
2. Press and hold the CAL key until displaying "CC out" and then release the key.
3. The balance displays "CC 0"
- If you want to change the calibration weight value, proceed to step 4.
- If you use the stored calibration weight value in the balance, proceed to step 5.
4. Press the RANGE key and adjust the calibration weight value using the following keys:
  - RE-ZERO key, the key to set the value of the digit selected.
  - RANGE key, the key to select the digit to change value.
  - PRINT key, the key to store a new weight value and return to step 3.
  - CAL key, the key to cancel this change and return to step 3.
5. Confirm that there is nothing on the pan and press the PRINT key. The balance measures the zero-point. Prevent vibration and drafts from affecting the balance.
6. Place the displayed target weight on the pan and press the PRINT key. The balance displays the measured target weight and displays it. Prevent vibration and drafts from affecting the balance.
7. Remove the weight from the pan after the balance displays "end".
8. If the "GLP output (info)" of the "Function Table " is set to "1" or "2", "Calibration Test Report" is output by the RS-232C interface.
9. The balance will automatically return to the weighing mode after calibration.








---

**ERROR CODES**

---

Display	Error Code	Description of the Error
	EC, E00	<b>Communications error</b> A protocol error occurred in communications. Confirm the format, baud rate and parity.
	EC, E01	<b>Undefined command error</b> An undefined command was received. Confirm command.
	EC, E02	<b>Not ready</b> A received command cannot be processed. Adjust the delay time to transmit the command. The balance received a <b>Q</b> command, but not in the weighing mode. The balance received a <b>Q</b> command while processing a RE-ZERO command.
	EC, E03	<b>Time over error</b> If t-Up 1 of the function list is set, the balance did not receive the next character of a command within the time limit of one second. Confirm communication.
	EC, E04	<b>Excess characters error</b> The balance received excessive characters in a command. Confirm command.
	EC, E06	<b>Format error</b> A command includes incorrect data. Confirm command. Data is numerically incorrect.
	EC, E07	<b>Range error for a parameter</b> The received data exceeds the range that the balance can accept. Confirm parameter range of command.
		<b>Internal condition information</b> There is no problem when this is displayed for a few seconds and returns to weighing mode. Turn the balance off then on again if this display is continuously displayed.
	EC, E11	<b>Stability error</b> The balance cannot stabilize due to an environmental problem. The balance is able to enter the weighing mode with the <b>CAL</b> command. Prevent vibration, drafts, temperature changes, static electricity and magnetic fields.
	EC, E16	<b>Internal weight error</b> This is a calibration error. Confirm that there is nothing on the pan and retry the calibration or calibration test.
	EC, E17	<b>Internal weight error</b> This is a calibration error. Retry the calibration or calibration test.

	EC, E20	<b>Calibration error</b> The calibration weight is too heavy. The balance is able to enter the weighing mode with the <b>CAL</b> command.
	EC, E21	<b>Calibration error</b> The calibration weight is too light. The balance is able to enter the weighing mode with the <b>CAL</b> command.
		<b>Overload</b> This is a warning that a weight beyond the balance capacity has been placed on the pan. Remove the weight from the pan.
		<b>Weighing pan error</b> This is a warning that the weight value is too light. Confirm that the weighing pan and the pan support are properly installed.
		<b>Unit weight, 100% weight error</b> The unit weight of the sample is very light in the counting mode, or the 100% sample is too light in percent mode. The balance cannot calculate it. Increase the unit weight or 100% weight.
		<b>ARA Zero error</b> The ARA (Automatic Response Adjustment) cannot be performed because there is something on the pan. Remove all matter from the pan. Press the <b>CAL</b> key to return to the weighing mode.
		<b>ARA Unstable error</b> The ARA (Automatic Response Adjustment) cannot be performed because of unstable weighing value. Correct the environment for the balance. Press the <b>CAL</b> key to return to the weighing mode.
		<b>Memory full</b> The memory data has reached 200 items. When data is deleted, new data can be stored.
		<b>Memory data error</b> The memory data is lost. Clear all memory data.
		<b>Unit weight information</b> This is advice regarding the sample number that is needed to set the unit weight. When the unit weight is computed and the sample number is too few, the required number is displayed for counting accuracy. Count and place the samples on the pan. Press the <b>PRINT</b> key to store the correct value.

## Legal-for-Trade Model Information

Model	GR-120-EC	GR-200-EC	GR-300-EC	GR-202-EC
Class	I	I	I	I
Max	120g	210g	310g	210g/42g
Min	10mg	10mg	10mg	10mg
Min Display d	0.1mg	0.1mg	0.1mg	0.1mg/0.01mg
Min Display e	1mg	1mg	1mg	1mg

Class	Item and Parameter	Description
dout Data output	ar-d Zero after output	Not Available Adjusts zero automatically after data is output
CS in Correction of internal weight		Not Available Correction of internal weight

- Calibration using an external weight (CalOut) is not available.
- Automatic span adjustment will be enabled (the balance will automatically calibrate using the internal weight due to environmental temperature changes).

**Note about Units:** When EU 2014/31/EU verification is performed only the following units will be available.

Unit Name	Unit Symbol	Unit Name	Unit Symbol
Grams / Milligram	g / mg	Pieces	pcs
Metric Caret	ct	Percent	%
Density	d		

### Note about Verification

#### Metrological note for EU countries

The identification plate of the scale has a CE mark followed by the metrology mark (M + year two-digit; framed) and the number of the notified body. Such a scale can be started up and operated on-site provided that it is not connected to an auxiliary device which is subject to metrological approval. Scales which are connected on site to an auxiliary device which is subject to metrological approval must undergo a conformity assessment or metrological approval. Only then, legal-for-trade use is permissible.

Scales without metrology mark or control marks may not be used in legal-for-trade applications.

Scales with an accuracy – Class II, with 'n' greater than 3300 have been adjusted to the place of use and may not be used elsewhere.

As per legal requirements, the scale user must ensure that the scale is used for its intended purpose. This includes the observance of legal requirements, especially in case of changes, amendments and extensions

#### Notes on verification

The verification of a scale is executed according to the respective statutory country regulations. The verification validity starts as soon as the device is put on the market.

---

### DISPOSAL OF THE DEVICE

---

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. To ensure that the device is disposed of correctly please recycle in accordance with local regulations. For further information contact the supplier from whom you purchased the product.

## EU Declaration of Conformity

This is a Declaration of Conformity for the instrument below.

INSTRUMENT MODEL	GR-EC Series
NAME AND ADDRESS OF THE MANUFACTURER	<b>A&amp;D Instruments Ltd</b> 24-26 Blacklands Way Abingdon Business Park Abingdon Oxfordshire OX14 1DY
TYPE APPROVAL CERTIFICATE	T5591
NOTIFIED BODY ISSUING THE APPROVAL FOR MODULE B: EU TYPE EXAMINATION (Annex II .1 of 2014/31)	NMI 0122
NOTIFIED BODY ISSUING THE APPROVAL FOR MODULE D: CONFORMITY TO TYPE BASED ON QUALITY ASSURANCE OF THE PRODUCTION PROCESS (Annex II.2 of 2014/31)	SGS UK 0120

Conformity to the following Directives is demonstrated by reference to the following harmonised standards or technical documentation:

APPLICABLE DIRECTIVES	HARMONISED STANDARDS OR OTHER TECHNICAL DOCUMENTATION
Directive on non-automatic weighing instrument (2014/31/EU)	EN45501:2015
Directive on Electromagnetic Compatibility (2014/30/EU)	EN 61326-1:2013 Electrical equipment for measurement control and laboratory use – EMC requirements Part 1: General Requirements
Directive on Electrical Equipment designed for use within certain voltage limits (2014/35/EU)	EN 60950-1:2006 +A11:2009+A1:2010+A12:2011+A2:2013 Safety of Information Technology Equipment
Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II) (2011/65/EU)	EN50581: 2012

This Declaration of Conformity is issued under the sole responsibility of the manufacturer:

Signed for and on behalf of:	<b>A&amp;D Instruments Ltd</b>
Place and date of issue:	Abingdon, UK, 1 <sup>st</sup> October 2018
Name, function:	J. Ghuman, Managing Director
Signature:	

---

**EU Declaration of Conformity**

---

**A&D Instruments Ltd. hereby declare that the following Weighing product conforms to the requirements of the council directives on ...**

Electromagnetic Compatibility (EMC) 2014/30/EU,  
Low Voltage Equipment (LVD) 2014/35/EU **and**  
Restriction of the use of certain Hazardous Substances (RoHS) 2011/65/EU

**provided that they bear the CE mark of conformity.**

**Model/Series....GR Series**

Standards applicable:

**EN 61326-1:2013**

Electrical equipment for measurement, control and laboratory use -EMC requirements  
Part 1: General requirements

**EN 60950-1:2006/A11:2009/A1:2010/A12:2011/A2:2013**

Safety of Information Technology Equipment

**EN 50581:2012**

Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

This Declaration of Conformity is issued under the sole responsibility of the manufacturer:  
A&D Instruments, 24-26 Blacklands Way, Abingdon, Oxfordshire, OX14 1DY, England

**CE Mark first applied 08 July 2002**

Signed for A&D Instruments on 1<sup>st</sup> October 2018

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